

Oversikt over akkrediterte analytter ved Avdeling for klinisk farmakologi sortert på analytt og matriks.

U= urin S= serum B= blod P=plasma

Analytt	Matriks	Analyseprinsipp
6-MAM (6-Monoacetylmorfin)	U, S, B	LC-MSMS
7-AF (7-aminoflunitrazepam)	U, B	LC-MSMS
7-AK (7-aminoklonazepam)	U, B	LC-MSMS
7-AN (7-aminonitrazepam)	U, B	LC-MSMS
Aceton	U, S, B	GC-FID
Alimemazin	S	LC-MSMS
Alprazolam	S, B	LC-MSMS
Amfetamin	U, S, B	LC-MSMS
Amisulprid	S	LC-MSMS
Amitriptylin	S	LC-MSMS
Apiksaban	P	LC-MSMS
Aripiprazol	S	LC-MSMS
Atomoksetin	S	LC-MSMS
Benzoyleggonin	U, S	LC-MSMS
Brivaracetam	S	LC-MSMS
Buprenorfin	U, S	LC-MSMS
Ciklosporin	B	LC-MSMS
Citalopram	S	SFC-MSMS
Dehydroaripiprazol	S	LC-MSMS
Deksamfetamin (D)	U, S	SFC-MSMS
Delorazepam	U	LC-MSMS
Desmetyldoksepin	S	LC-MSMS
Desmetylklomipramin	S	LC-MSMS
Diazepam	S, B	LC-MSMS
Digitoksin	S	CEDIA, AU680 (autoanalyzer)
DMD (desmetyldiazepam)	U, S, B	LC-MSMS
DM-venlafaksin	S	LC-MSMS
Dokseptin	S	LC-MSMS
Duloksetin	S	LC-MSMS
Escitalopram	S	SFC-MSMS
Edoksaban	P	LC-MSMS
Etanol	U, S, B	GC-FID
Etosuksimid	S	LC-MSMS
Etylenglykol	S, B	GC-MS
Etylglukuronid	U	LC-MSMS
Etylglukuronid	U	DRI, AU680 (autoanalyzer)
Etylmorfin	U, S, B	LC-MSMS
Etylsulfat	U	LC-MSMS

Everolimus	B	LC-MSMS
Fenobarbital	S	CEDIA, AU680 (autoanalyzer)
Fentanyl	U, S	LC-MSMS
Fenytoin	S	CEDIA, AU680 (autoanalyzer)
Flekainid	S	LC-MSMS
Flunitrazepam	S, B	LC-MSMS
Fluoksetin	S	LC-MSMS
Flupentisol	S	LC-MSMS
Fluvoksamin	S	LC-MSMS
Fosfatidyletanol (PEth)	B	LC-MSMS
Gabapentin	S	LC-MSMS
GHB (Gammahydroksybutyrat)	U	LC-MSMS
Haloperidol	S	LC-MSMS
Isopropanol	U, S, B	GC-FID
Karbamazepin	S	LC-MSMS
Karbamazepinepoksid	S	LC-MSMS
Ketamin	U	LC-MSMS
Ketobemidon	U, S	LC-MSMS
Klomipramin	S	LC-MSMS
Klonazepam	S, B	LC-MSMS
Klorprotiksen	S	LC-MSMS
Klozapin	S	LC-MSMS
Kodein	U, S, B	LC-MSMS
Kodein-6-glukuronid (K6G)	U, S, B	LC-MSMS
Koffein	S	LC-MSMS
Kreatinin	U	DRI, AU680 (autoanalyzer)
Lakosamid	S	LC-MSMS
Lamotrigin	S, B	LC-MSMS
Levetiracetam	S	LC-MSMS
Levoamfetamin (L)	U, S	SFC-MSMS
Levomepromazin	S	LC-MSMS
Likarbazepin	S	LC-MSMS
Litium	S	Fotometri, AU680 (autoanalyzer)
Lorazepam	U, S, B	LC-MSMS
Morfin-3-glukuronid (M3G)	U, S, B	LC-MSMS
Morfin-6-glukuronid (M6G)	U, S, B	LC-MSMS
MDA (4-methyl-dioxy-amfetamin)	U, S, B	LC-MSMS
MDMA (3,4-Metylendioksymetamfetamin)	U, S, B	LC-MSMS
Metadon	U, S, B	LC-MSMS
Metamfetamin	U, S, B	LC-MSMS
Metanol	U, S, B	GC-FID
Mianserin	S	LC-MSMS
Midazolam	S, B	LC-MSMS
Mirtazapin	S	LC-MSMS

Moklobemid	S	LC-MSMS
Morfin	U, S, B	LC-MSMS
Mykofenolsyre	S, P	LC-MSMS
Nitrazepam	S, B	LC-MSMS
Norfluoksetin	S	LC-MSMS
Nortriptylin	S	LC-MSMS
O-desmetyltramadol	U	LC-MSMS
OH-alprazolam	U	LC-MSMS
OH-bupropion	S	LC-MSMS
Oksazepam	U, S, B	LC-MSMS
Oksykodon	U, S, B	LC-MSMS
Olanzapin	S	LC-MSMS
Paliperidon (OH-risperidon)	S	LC-MSMS
Paracetamol	S	DRI, AU680 (autoanalyzer)
Paroksetin	S	LC-MSMS
Perampanel	S	LC-MSMS
Perfenazin	S	LC-MSMS
Petidin	U	LC-MSMS
PMMA (Parametoksymetamfetamin)	U	LC-MSMS
Pregabalin	U, S	LC-MSMS
Quetiapin	S	LC-MSMS
Risperidon	S	LC-MSMS
Ritalinsyre	U, S, B	LC-MSMS
Rivaroksaban	P	LC-MSMS
Rufinamid	S	LC-MSMS
Sertindol	S	LC-MSMS
Sertralin	S	LC-MSMS
Takrolimus	B	LC-MSMS
Tapentadol	U	LC-MSMS
Temazepam	U	LC-MSMS
Teofyllin	S	LC-MSMS
THCA	U	LC-MSMS
Topiramat	S	LC-MSMS
Tramadol	U, S	LC-MSMS
Trimipramin	S	LC-MSMS
Valproat	S	CEDIA, AU680 (autoanalyzer)
Venlafaksin	S	LC-MSMS
Vortioksetin	S	LC-MSMS
Ziprasidon	S	LC-MSMS
Zolpidem	U, S	LC-MSMS
Zonisamid	S	LC-MSMS
Zopiklon	U, S, B	LC-MSMS
Zuklopentiksol	S	LC-MSMS

