

Avdeling for medisinsk mikrobiologi (AMM) – Infeksjonsimmunologi (serologi)			
Cytomegalovirus	Parotittvirus (kusma)	Hepatitt A-virus	Yersinia enterocolitica
Epstein-Barr-virus	Morbillivirus (meslinger)	Hepatitt B-virus (akt. inf./bærer)	Streptokokk-anti-DNase B
Herpes simplex-virus	Rubellavirus	Anti-HBs (vaksinetiter)	AST
Varizella-zoster-virus	Parvovirus B19		Helicobacter pylori
Puumala/Hantavirus	Toxoplasma gondii	Hepatitt C-virus	Francisella tularensis
Denguevirus		HIV	Borrelia burgdorferi
		Syfilis	Bordetella pertussis (kikhoste)



Hver fargegruppe trenger alle hvert sitt **serum-gelrør** 5 mL.



UNNTAK

*PCR- us:


Obs!
Begrenset holdbarhet

*Cytomegalovirus, Epstein-Barr-virus, Humant herpes virus 6.

*Hiv, Hepatitt B og C, Parvovirus B19



EDTA-rør



EDTA-rør m/gel (PPT-rør)

****QuantiFERON:** Spesialglass, oppbevares i romtemperatur

Se elektronisk brukerhåndbok for spesiell prøvebehandling.

www.stolav.no/lab - Avdeling for medisinsk mikrobiologi)

Ved spørsmål kontakt vakthavende mikrobiolog tlf. 95 97 94 97